# DR. MÖLLER & SCHMELZ GmbH

# Gesellschaft für angewandte Mikrobiologie



# **VRBD-NPS**

Version: 10/2017

M&S item numbers: 1225 (50 / PK) and 1225-H (100 / PK)

Profile: Dehydrated nutrient pad sets 50 mm in petri dishes, sterile

Color: Light red

Storage: Dark and dry at room temperature

Shelf life: 2 years after sterilization

## Description and application range

VRBD-NPS are used for the determination of enteric bacteria from water, food and other samples. One characteristic for enteric bacteria is the formation of acids from dextrose, which is indicated by red colonies due to the neutral red indicator. Growth of accompanying microorganisms is inhibited by crystal violet and bile salts. Note: Growth on this medium can only be taken as hint. For confirmation further tests need to be performed. The medium is manufactured and quality tested in compliance with ISO 11133:2014 standard.

### **Typical composition**

Enzymatic digest of animal tissue	7.0 g/l
Yeast extract	3.0 g/l
Sodium chloride	5.0 g/l
Bile salts	1.5 g/l
Dextrose	10.0 g/l
Neutral red	0.05 g/l
Crystal violet	0.002 g/l

Final pH: 7.2 ± 0.2 at 25 °C

#### Microbiological quality control

#### **Bacterial contamination**

Incubation: aerobically at room temperature for 3 days, specification: no growth

#### **Productivity** quantitative analysis

Incubation: aerobically at 37  $\pm$  1 °C for 24  $\pm$  2 h, approx. inoculum: 80 – 120 CFU

Microorganism	Test strain	Specification	Appearance
Escherichia coli	WDCM 00012	P <sub>R</sub> ≥ 0,5	Red
Enterobacter aerogenes	WDCM 00175	P <sub>R</sub> ≥ 0,5	Light red